

Development of Laser Induced Breakdown Spectroscopy Sensor for Molten Salt Reactor Off-Gas Stream


Hunter Andrews

K. Myhre, Joanna McFarlane, Shay Chapel, N. Dianne Ezell

MSR Workshop 2021

October 13th, 2021


We're seeking to enable molten salt reactor (MSR) deployment by developing technology for off-gas analysis



- Quantify fission and activation products in off-gas system
- Monitor off-gas treatment component efficiency

What is laser-induced breakdown spectroscopy (LIBS)?


Laser ablation and optical emission collection



Why laser induced breakdown spectroscopy (LIBS)?

Elemental analysis via optical emission spectroscopy of laser induced plasmas

- Benefits:
 - Little to no sample prep
 - Useful for solid, liquid, and gas analysis
 - Typical sensitivities of ppm
 - Quasi-nondestructive (nanograms per shot)
 - Can be completely fiber optic based





Design Challenges

- Producing a surrogate aerosol stream
- Designing a measurement cell








A collision nebulizer was selected to generate aerosol stream for analysis




[1] A. Williams, Virginia Commonwealth University, 2016.


Aerosol particle sizes range from ~1-10 μm in diameter




A sheath gas approach was used to contain aerosol stream during measurement




Sheath gas system successfully contained aerosol stream



Sheath gas is turned on and off repeatedly



A complete prototype aerosol measurement system was manufactured





Aqueous aerosol monitoring system

- The three representative elements selected were Gd, Nd, and Sm
 - All act as neutron poisons in a reactor making their concentrations of interest to operators.
 - Concentrations ranging from 0 to 2000 ppm in liquid reservoir
 - 2000 ppm in reservoir \approx 5 ppb in aerosol stream



Regression models were developed using a set of calibration samples




Sample 2 is Gd dominant
Sample 11 is Nd dominant
Sample 13 is Sm dominant

What is PLS regression?


Partial Least Squares Regression (PLS)

1. Matrices X and Y are decomposed into latent structures in an iterative process.
2. The latent structures corresponding to the most variation of Y (u_i) is explained by a latent structure in X (t_i) which explains it the best.




Note: the goal is to explain the most variance in Y , not necessarily X

Predicted concentrations of validation samples match ICP-OES measurements




Performing a demonstration to show the measurement system's capability to provide real-time monitoring

- Stock solutions pumped in/out of reservoir to nebulizer to allow real-time changes in concentrations




Initial concentration predictions leave room for improvement






The RMSE values were calculated to be 249, 108, and 99.8 ppm for Gd, Nd, and Sm, respectively.

A genetic algorithm is an optimization approach based on Darwin's theory of evolution



A genetic algorithm was used to refine PLS feature selection




Number of LVs

Gd: 5 → 4

Nd: 6 → 3

Sm: 9 → 3

Genetic algorithm filtered models showed dramatically improved results



This corresponds to a 73, 18, and 25% decrease in RMSE value from the previous model predictions.



The RMSE values were calculated to be 66.5, 89.2, and 75.3 ppm for Gd, Nd, and Sm, respectively.

Gaseous monitoring system

- The four representative elements selected were Xe, Kr, Cs, and Rb
 - Xe and Kr are expected to be large contributors to the off-gas loading
 - 0 to 10 g/h
 - Cs and Rb are the corresponding daughters of these gases
 - 0 to 2000 ppm




Krypton and xenon showed strong spectral responses





Collected LIBS spectrum of sample 7 containing 1.21 and 1.90 mol% of Xe and Kr and 2000 ppm Rb.

A closer look at collected spectra reveal strong gas peaks




Collected LIBS spectrum of a blank sample (sample 5) and sample 14 containing 1.21 and 1.90 mol% of Xe and Kr and 2000 and 1800.9 ppm of Cs and Rb.

Xe and Kr univariate models were developed



Univariate calibration models, dashed lines, comparing (a) 881.94 nm Xe I peak area and (b) 760.15 nm Kr I peak area response to changes to Xe and Kr mass flow rates, respectively.




Analytes in test samples were predicted using constructed models



*signifies a univariate regression model

The Kr concentration was successfully predicted in real-time

Model predictions for the concentration of Kr gas compared to converted flow meter measurements in systems where the gas flow rate is changed




In summary, we have demonstrated a LIBS capability for monitoring aerosols and gases.

- Our sheathed gas measuring system allows effective measurement of the sample stream while protecting optical components
- Multiple elements were simultaneously quantified in the two studies discussed
- A genetic algorithm was used to refine a quantitative model
- Both gases and aerosols were able to be monitored in real-time

Where is this project headed next?

- Future work involves completion of a molten salt aerosol test stand and integrating LIBS sensors with off-gas treatment systems.



Acknowledgments

MSR Campaign Team

NTD for Technology Development

NTD for MSR Licensing

Licensing Framework

Technology Development

Off-gas System

Chemistry

Tritium

Thermophysical Properties

Structural Materials

Graphite

Advanced Material Development

Waste forms

Salt Spill

Patricia Paviet

Jim Kinsey

David Holcomb, Dave Luxat

Kevin Robb, Jordan Massengale

Joanna McFarlane, Hunter Andrews, Shay Chapel

Amanda Lines, Sam Bryan, Heather Felmy, Adan Medina

Jake McMurray, Mark Williamson, Toni Karlsson, Nathan Hoyt

Paul Humrickhouse

Melissa Rose, Ryan Gallagher, Dianne Ezell, Marissa Monreal

Bruce Pint, Jim Keiser

Nidia Gallego

Sam Sham

Brian Riley

Sara Thomas, Bill Ebert

Questions?

E-mail: andrewshb@ornl.gov



End